


easyvision Soria Multifocal Fitting Guide



easyvision Soria Multifocal has high 1st time fitting success of 95% and 100% within 2 lenses¹

Neuro Adaptive™ Design Centre near vision



Easy adaption with natural smooth transitions through near, intermediate and far viewing.

Successful fitting with easyvision Soria Multifocal

easyvision Soria Multifocal has been specifically designed to make Multifocal fitting easier.

- Maximise the plus sphere powers during distance vision refraction
- Account for vertex distance (>+/-4D) and mean spherical equivalent for astigmatism up to 1D
- Determine dominant eye
- Select the initial add. Extra distance plus power should be offered using real vision assessments as first choice and is often well accepted.

Initial lens selection	Distance vision enhancements	Near vision enhancements
Add +0.75D to +1.25D 	Refine distance power 	Increase plus power
Add +1.50D to +1.75D 	Refine distance power 	Switch to HIGH
Add +2.00D and above 	Refine distance power 	Increase plus power

General tip:
Always push the plus giving the least reading Add needed. For enhancements use full aperture trial lenses increasing the power as needed. Binocular adjustments are also possible.

*Remove extra +0.50D from initial lens selection

¹Menicon data on file 2021 from an Independent clinical study.